

PRB-HLD

Track Stand Probe Holder



The **PRB-HLD** mounts the fiber optic probe that is included as a standard design feature on the PeakSeeker™ and PeakSeeker Pro™ Raman Spectrometer Systems.

The PRB-HLD is used for Raman analysis of samples placed on its base plate. A rack and pinion track stand provides fine vertical adjustment for precision Raman probe positioning.

Specifications

Weight	6.6 lbs (3.0 kg)
Size of the base	12.6 × 12.0 × 0.63 inches (320 × 305 × 16 mm)
Travel distance	250 mm
Vertical length	325 mm

*Specifications subject to change without notice.
©2010 All rights reserved.*

